

# Human Symmetric dimethylarginine (SMDA) ELISA Kit

Catalog No: #EK20651

Package Size: #EK20651-1 48T #EK20651-2 96T

## Description

|                    |                                                   |
|--------------------|---------------------------------------------------|
| Product Name       | Human Symmetric dimethylarginine (SMDA) ELISA Kit |
| Applications       | ELISA                                             |
| Species Reactivity | Human                                             |

## Application Details

Detect Range:3.12-200 ng/mL

Sensitivity:1.19ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.