

Fish Thyroid stimulating hormone (TSH) ELISA Kit

Catalog No: #EK20652

Package Size: #EK20652-1 48T #EK20652-2 96T



Description

Product Name	Fish Thyroid stimulating hormone (TSH) ELISA Kit
Applications	ELISA
Species Reactivity	Fish

Application Details

Detect Range:0.3-12 μ g/mL

Sensitivity:0.1 μ g/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.