

# Human Ig alpha-1 chain C region (IGHA1) ELISA Kit

Catalog No: #EK20653

Package Size: #EK20653-1 48T #EK20653-2 96T



## Description

|                    |                                                   |
|--------------------|---------------------------------------------------|
| Product Name       | Human Ig alpha-1 chain C region (IGHA1) ELISA Kit |
| Applications       | ELISA                                             |
| Species Reactivity | Human                                             |

## Application Details

Detect Range:1.56-100&micro;g/mL

Sensitivity:0.52&micro;g/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.