

Rat SH3 and multiple ankyrin repeat domains protein 3 (SHANK3) ELISA Kit

Catalog No: #EK20661

Package Size: #EK20661-1 48T #EK20661-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat SH3 and multiple ankyrin repeat domains protein 3 (SHANK3) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:23.4pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.