

Human Protein S100-A3 (S100A3) ELISA Kit

Catalog No: #EK20663

Package Size: #EK20663-1 48T #EK20663-2 96T



Description

Product Name	Human Protein S100-A3 (S100A3) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:1.56-100ng/mL

Sensitivity:0.57ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.