

Monkey Renin (REN) ELISA Kit

Catalog No: #EK20667

Package Size: #EK20667-1 48T #EK20667-2 96T



Description

Product Name	Monkey Renin (REN) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:31.2-2000 pg/mL

Sensitivity:12.5pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.