

Rat Actin, aortic smooth muscle (ACTA2) ELISA Kit

Catalog No: #EK20673

Package Size: #EK20673-1 48T #EK20673-2 96T

Description

Product Name	Rat Actin, aortic smooth muscle (ACTA2) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:6.25-400 ng/mL

Sensitivity:2.87ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.