

Mouse Lipid phosphate phosphohydrolase 3 (PLPP3) ELISA Kit

Catalog No: #EK20678

Package Size: #EK20678-1 48T #EK20678-2 96T

Description

Product Name	Mouse Lipid phosphate phosphohydrolase 3 (PLPP3) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000 pg/mL

Sensitivity:10.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.