

# Mouse Dopa Decarboxylase (DDC) ELISA Kit

Catalog No: #EK20679

Package Size: #EK20679-1 48T #EK20679-2 96T

## Description

|                    |                                          |
|--------------------|------------------------------------------|
| Product Name       | Mouse Dopa Decarboxylase (DDC) ELISA Kit |
| Applications       | ELISA                                    |
| Species Reactivity | Mouse                                    |

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.062ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.