

Human Quinone oxidoreductase PIG3 (TP53I3) ELISA Kit

Catalog No: #EK20681

Package Size: #EK20681-1 48T #EK20681-2 96T

Description

Product Name	Human Quinone oxidoreductase PIG3 (TP53I3) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:34pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.