

Human Leucine-rich repeat transmembrane protein FLRT2 (FLRT2) ELISA Kit

Catalog No: #EK20683

Package Size: #EK20683-1 48T #EK20683-2 96T

Description

Product Name	Human Leucine-rich repeat transmembrane protein FLRT2 (FLRT2) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:32.6pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture
supernates or other biological fluids.Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.