

Monkey Hemoglobin subunit Beta (HBB) ELISA Kit

Catalog No: #EK20686

Package Size: #EK20686-1 48T #EK20686-2 96T

Description

Product Name	Monkey Hemoglobin subunit Beta (HBB) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:0.78-50 ng/mL

Sensitivity:0.31ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.