

Human 45 kDa calcium-binding protein (SDF4) ELISA Kit

Catalog No: #EK20690

Package Size: #EK20690-1 48T #EK20690-2 96T

Description

Product Name	Human 45 kDa calcium-binding protein (SDF4) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.067ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.