

# Monkey Apolipoprotein C3 (APOC3) ELISA Kit

Catalog No: #EK20694

Package Size: #EK20694-1 48T #EK20694-2 96T



## Description

|                    |                                            |
|--------------------|--------------------------------------------|
| Product Name       | Monkey Apolipoprotein C3 (APOC3) ELISA Kit |
| Applications       | ELISA                                      |
| Species Reactivity | Monkey                                     |

## Application Details

Detect Range:1.56-100 ng/mL

Sensitivity:0.67ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.