

Canine C-Reactive Protein (CRP) ELISA Kit

Catalog No: #EK20697

Package Size: #EK20697 48T

Description

Product Name	Canine C-Reactive Protein (CRP) ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range:0.312-20 μ g/mL

Sensitivity:0.105 μ g/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.