

Rat Immunoglobulin G4 (IgG4) ELISA Kit

Catalog No: #EK20700

Package Size: #EK20700-1 48T #EK20700-2 96T



Description

Product Name	Rat Immunoglobulin G4 (IgG4) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:15.6-1000 ng/mL

Sensitivity:6.7ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.