

Rat Anti-cyclic Citrullinated Peptide Antibody (CCP-Ab) ELISA Kit

Catalog No: #EK20703

Package Size: #EK20703-1 48T #EK20703-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Anti-cyclic Citrullinated Peptide Antibody (CCP-Ab) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:3.12-200 IU/mL

Sensitivity:1.27 IU/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.