

# Canine Insulin (INS) ELISA Kit

Catalog No: #EK20712

Package Size: #EK20712-1 48T #EK20712-2 96T

## Description

|                    |                                |
|--------------------|--------------------------------|
| Product Name       | Canine Insulin (INS) ELISA Kit |
| Applications       | ELISA                          |
| Species Reactivity | Canine                         |

## Application Details

proformaDetect Range:3.12-200ng/mL

Sensitivity:1.16ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.