

Porcine Secretory immunoglobulin A (SIgA) ELISA Kit

Catalog No: #EK20715

Package Size: #EK20715-1 48T #EK20715-2 96T

Description

Product Name	Porcine Secretory immunoglobulin A (SIgA) ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

Application Details

Detect Range:1.56-100 μ g/mL

Sensitivity:0.51 μ g/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.