

Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit

Catalog No: #EK20729

Package Size: #EK20729-1 48T #EK20729-2 96T

Description

Product Name	Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:12.5-800pg/mL

Sensitivity:5.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.