

Mouse Chemokine C-X-C-Motif Receptor 6 (CXCR6)
ELISA Kit

Catalog No: #EK20730

Package Size: #EK20730-1 48T #EK20730-2 96T



Description

Product Name	Mouse Chemokine C-X-C-Motif Receptor 6 (CXCR6) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.058ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.