

Monkey Complement Component 4 (C4) ELISA Kit

Catalog No: #EK20736

Package Size: #EK20736-1 48T #EK20736-2 96T

Description

Product Name	Monkey Complement Component 4 (C4) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:62.5-4000ng/mL

Sensitivity:19.3ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.