

Human Immunoglobulin Lambda Like Polypeptide 5 (IgIL5) ELISA Kit

Catalog No: #EK20737

Package Size: #EK20737-1 48T #EK20737-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Immunoglobulin Lambda Like Polypeptide 5 (IgIL5) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL

Sensitivity:0.092ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.