

Human Elastase (ELA) ELISA Kit

Catalog No: #EK20741

Package Size: #EK20741-1 48T #EK20741-2 96T

Description

Product Name	Human Elastase (ELA) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL

Sensitivity:0.107ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.