

Monkey Amyloid Beta Peptide 1-40 (Ab1-40) ELISA Kit

Catalog No: #EK20747

Package Size: #EK20747-1 48T #EK20747-2 96T

Description

Product Name	Monkey Amyloid Beta Peptide 1-40 (Ab1-40) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:125-8000pg/mL

Sensitivity:43pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.