

Human Calcineurin-like phosphoesterase domain-containing protein 1 (CPPED1) ELISA Kit

Catalog No: #EK20748

Package Size: #EK20748-1 48T #EK20748-2 96T

Description

Product Name	Human Calcineurin-like phosphoesterase domain-containing protein 1 (CPPED1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.064ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.