

# Human Coiled-coil domain-containing protein 134 (CCDC134) ELISA Kit

Catalog No: #EK20749

Package Size: #EK20749-1 48T #EK20749-2 96T

## Description

|                    |                                                                     |
|--------------------|---------------------------------------------------------------------|
| Product Name       | Human Coiled-coil domain-containing protein 134 (CCDC134) ELISA Kit |
| Applications       | ELISA                                                               |
| Species Reactivity | Human                                                               |

## Application Details

|                                                                                                    |
|----------------------------------------------------------------------------------------------------|
| Detect Range:0.156-10ng/mL                                                                         |
| Sensitivity:0.067ng/mL                                                                             |
| Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids |
| Sample Volume: 50 &mu;L                                                                            |
| Assay Time:1-4.5h                                                                                  |
| Detection wavelength:450 nm                                                                        |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.