

Human N-acetylglucosamine-1-phosphotransferase subunit gamma (GNPTG) ELISA Kit

Catalog No: #EK20759

Package Size: #EK20759-1 48T #EK20759-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human N-acetylglucosamine-1-phosphotransferase subunit gamma (GNPTG) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL

Sensitivity:0.103ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.