

Monkey Neurofilament protein L (NF-L) ELISA Kit

Catalog No: #EK20762

Package Size: #EK20762-1 48T #EK20762-2 96T

Description

Product Name	Monkey Neurofilament protein L (NF-L) ELISA Kit
Applications	ELISA
Species Reactivity	Monkey

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.062ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.