

Human Laminin Subunit Alpha-1 (LAMA1) ELISA Kit

Catalog No: #EK20764

Package Size: #EK20764-1 48T #EK20764-2 96T



Description

Product Name	Human Laminin Subunit Alpha-1 (LAMA1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL
Sensitivity:31pg/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.