

Human Interferon-inducible double stranded
RNA-dependent protein kinase activator A (PRKRA)
ELISA Kit

Catalog No: #EK20766

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #EK20766-1 48T #EK20766-2 96T

Description

Product Name	Human Interferon-inducible double stranded RNA-dependent protein kinase activator A (PRKRA) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:125-8000pg/ml
Sensitivity:43pg/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.