

# Human Haloacid dehalogenase-like hydrolase domain-containing protein 2 (HDHD2) ELISA Kit

Catalog No: #EK20768

Package Size: #EK20768-1 48T #EK20768-2 96T

## Description

|                    |                                                                                          |
|--------------------|------------------------------------------------------------------------------------------|
| Product Name       | Human Haloacid dehalogenase-like hydrolase domain-containing protein 2 (HDHD2) ELISA Kit |
| Applications       | ELISA                                                                                    |
| Species Reactivity | Human                                                                                    |

## Application Details

Detect Range:125-8000pg/mL

Sensitivity:49pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.