

# Monkey SPARC-related modular calcium-binding protein 1 (SMOC1) ELISA Kit

Catalog No: #EK20769

Package Size: #EK20769-1 48T #EK20769-2 96T

## Description

|                    |                                                                          |
|--------------------|--------------------------------------------------------------------------|
| Product Name       | Monkey SPARC-related modular calcium-binding protein 1 (SMOC1) ELISA Kit |
| Applications       | ELISA                                                                    |
| Species Reactivity | Monkey                                                                   |

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.064ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.