

Mouse UTP--glucose-1-phosphate uridylyltransferase (UGP2) ELISA Kit

Catalog No: #EK20781

Package Size: #EK20781-1 48T #EK20781-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse UTP--glucose-1-phosphate uridylyltransferase (UGP2) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.7pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.