

Bovine Immunoglobulin A (IgA) ELISA Kit

Catalog No: #EK20787

Package Size: #EK20787-1 48T #EK20787-2 96T



Description

Product Name	Bovine Immunoglobulin A (IgA) ELISA Kit
Applications	ELISA
Species Reactivity	Bovine

Application Details

Detect Range:31.2-2000μg/mL

Sensitivity:12.6μg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.