

Porcine Beta-2-Microglobulin (b2M) ELISA Kit

Catalog No: #EK20789

Package Size: #EK20789-1 48T #EK20789-2 96T

Description

Product Name	Porcine Beta-2-Microglobulin (b2M) ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.27ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.