

Mouse Caspase 11 (CASP11) ELISA Kit

Catalog No: #EK20799

Package Size: #EK20799-1 48T #EK20799-2 96T

Description

Product Name	Mouse Caspase 11 (CASP11) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.059ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.