

Rat Plasmin-antiplasmin complex (PAP) ELISA Kit

Catalog No: #EK20805

Package Size: #EK20805-1 96T

Description

Product Name	Rat Plasmin-antiplasmin complex (PAP) ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:3.12-200ng/mL

Sensitivity:1.56ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.