

Human 39S ribosomal protein L21, mitochondrial (MRPL21) ELISA Kit

Catalog No: #EK20809

Package Size: #EK20809-1 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human 39S ribosomal protein L21, mitochondrial (MRPL21) ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	L21mt; MGC62013; MRP-L21; 39S ribosomal protein L21; mitochondrial

Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.