

# Human 28S ribosomal protein S29, mitochondrial (DAP3) ELISA Kit

Catalog No: #EK20821

Package Size: #EK20821-1 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                                                                           |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Human 28S ribosomal protein S29, mitochondrial (DAP3) ELISA Kit                                                                           |
| Applications       | ELISA                                                                                                                                     |
| Species Reactivity | Human                                                                                                                                     |
| Other Names        | DAP-3; DKFZp686G12159; MGC126058; MGC126059; MRP-S29; MRPS29; bMRP-10; death-associated protein 3 mitochondrial 28S ribosomal protein S29 |

## Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.