

Mouse Tetraspanin-4(TSPAN4) ELISA Kit

Catalog No: #EK20828

Package Size: #EK20828-1 48T #EK20828-2 96T

Description

Product Name	Mouse Tetraspanin-4(TSPAN4) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/ml

Sensitivity:12.5pg/ml

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.