

Mouse Forkhead box protein P1(FOXP1) ELISA Kit

Catalog No: #EK20830

Package Size: #EK20830-1 48T #EK20830-2 96T



Description

Product Name	Mouse Forkhead box protein P1(FOXP1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.062ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.