

Human Nucleotide triphosphate diphosphatase NUDT15(NUDT15) ELISA Kit

Catalog No: #EK20839

Package Size: #EK20839-1 48T #EK20839-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Nucleotide triphosphate diphosphatase NUDT15(NUDT15) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.064ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.