

Bovine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit

Catalog No: #EK20841

Package Size: #EK20841-1 48T #EK20841-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit
Applications	ELISA
Species Reactivity	Bovine

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:18.8pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.