

Human TRAF-type zinc finger domain-containing protein 1(TRAFFD1) ELISA Kit

Catalog No: #EK20842

Package Size: #EK20842-1 48T #EK20842-2 96T



Description

Product Name	Human TRAF-type zinc finger domain-containing protein 1(TRAFFD1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.069ng/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.