

Human Large ribosomal subunit protein bL27m(MRPL27) ELISA Kit

Catalog No: #EK20845

Package Size: #EK20845-1 48T #EK20845-2 96T

Description

Product Name	Human Large ribosomal subunit protein bL27m(MRPL27) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.35ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.