

# Mouse Myelin Oligodendrocyte Glycoprotein (MOG) ELISA Kit

Catalog No: #EK20850

Package Size: #EK20850-1 48T #EK20850-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                           |
|--------------------|-----------------------------------------------------------|
| Product Name       | Mouse Myelin Oligodendrocyte Glycoprotein (MOG) ELISA Kit |
| Applications       | ELISA                                                     |
| Species Reactivity | Mouse                                                     |

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:10.7pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &amp;mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.