

# Human Cadherin-related family member 2 (CDHR2) ELISA Kit

Catalog No: #EK20852

Package Size: #EK20852-1 48T #EK20852-2 96T

## Description

|                    |                                                          |
|--------------------|----------------------------------------------------------|
| Product Name       | Human Cadherin-related family member 2 (CDHR2) ELISA Kit |
| Applications       | ELISA                                                    |
| Species Reactivity | Human                                                    |

## Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.062ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.