

Mouse Muscle-specific RING finger protein 1(Murf1) ELISA Kit

Catalog No: #EK20854

Package Size: #EK20854-1 48T #EK20854-2 96T

Description

Product Name	Mouse Muscle-specific RING finger protein 1(Murf1) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.054ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

References

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.