

Human Leukocyte immunoglobulin-like receptor subfamily B member 1(LILRB1) ELISA Kit

Catalog No: #EK20863

Package Size: #EK20863-1 48T #EK20863-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Leukocyte immunoglobulin-like receptor subfamily B member 1(LILRB1) ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:4.75pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.